


Jeio Tech Stainless Steel Flask Clamp (2 Litres)

Product Code: 23555

Brand: Jeio Tech

The Jeio Tech Stainless Steel Flask Clamp (2 Litres) is designed for 2 litre Erlenmeyer flasks when using laboratory incubators and shakers.



Specifications:

Erlenmeyer Flask Size		2
Max. Mountable Flask Clamps	IST-3075(R)	4
	IST-4075(R)	5
	ISF-7100(R)	11
	ISF-7200(R)	
	ISS-7100(R)	
	ISS-7200(R)	
	ISS-3075(R)	4
	ISS-4075(R)	6
	OS-2000	-
	OS-3000	4

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com