


## Jeio Tech Stainless Steel Flask Clamp (300ml)

**Product Code:** 23556

**Brand:** Jeio Tech

The Jeio Tech Stainless Steel Flask Clamp (300ml) is designed for 300ml Erlenmeyer flasks when using laboratory incubators and shakers.



### Specifications:

Erlenmeyer Flask Size (ml)		300
<b>Max. Mountable Flask Clamps</b>	<b>IST-3075(R)</b>	14
	<b>IST-4075(R)</b>	20
	<b>ISF-7100(R)</b>	39
	<b>ISF-7200(R)</b>	
	<b>ISS-7100(R)</b>	
	<b>ISS-7200(R)</b>	
	<b>ISS-3075(R)</b>	14
	<b>ISS-4075(R)</b>	24
	<b>OS-2000</b>	8
	<b>OS-3000</b>	14

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)