

Jeio Tech Stainless Steel Flask Clamp (500ml)

Product Code: 23553

Brand: Jeio Tech

The Jeio Tech Stainless Steel Flask Clamp (500ml) is designed for 500ml Erlenmeyer flasks when using laboratory incubators and shakers.



Specifications:

Erlenmeyer Flask Size (ml)		500
Max. Mountable Flask Clamps	IST-3075(R)	9
	IST-4075(R)	16
	ISF-7100(R)	28
	ISF-7200(R)	
	ISS-7100(R)	
	ISS-7200(R)	
	ISS-3075(R)	9
	ISS-4075(R)	16
	OS-2000	6
	OS-3000	10

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com