


## Jeio Tech Stainless Steel Flask Clamp (6 Litres)

**Product Code:** 23559

**Brand:** Jeio Tech

The Jeio Tech Stainless Steel Flask Clamp (6 Litres) is designed for 6 litre Erlenmeyer flasks when using laboratory incubators and shakers.



### Specifications:

<b>Erlenmeyer Flask Size</b>		6
<b>Max. Mountable Flask Clamps</b>	<b>IST-3075(R)</b>	-
	<b>IST-4075(R)</b>	-
	<b>ISF-7100(R)</b>	5
	<b>ISF-7200(R)</b>	
	<b>ISS-7100(R)</b>	-
	<b>ISS-7200(R)</b>	
	<b>ISS-3075(R)</b>	-
	<b>ISS-4075(R)</b>	2
	<b>OS-2000</b>	-
	<b>OS-3000</b>	1

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)