

Jeio Tech Test Tube Rack for 10 Test Tubes (25mm)

Product Code: AAA41539

Brand: Jeio Tech

Jeio Tech's Test Tube Rack for 10 Test Tubes (25mm) is highly durable, and features adjustable angle and tube height support without requiring any base plate. It can be easily mounted and released.

Suitable for Jeio Tech 3.5 Litre Water Bath (BW-05H) and Jeio Tech 3.5 Litre Water Bath (BW-05B).



Specifications:

Diameter (mm)	25
Insert Capacity of Racks	BW-05H BW-05B
	1 ea

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com