

## Jeio Tech Test Tube Rack for 176 Test Tubes (10mm)

**Product Code:** AAA45552

**Brand:** Jeio Tech

Jeio Tech's Test Tube Rack for 176 Test Tubes (10mm) is highly durable and features adjustable angle and tube height support without requiring any base plate. It can easily be mounted and released.

Suitable for Jeio Tech Water Baths (BW-H Series) and Jeio Tech Water Baths (BW-B Series).



### Specifications:

<b>Diameter (mm)</b>	10	
<b>Insert Capacity of Racks</b>	BW-10H BW-10B	BW-20H BW-20B
	1ea	2ea

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)