

Jeio Tech Test Tube Rack for 32 Test Tubes (16mm)

Product Code: AAA44585

Brand: Jeio Tech

Jeio Tech's Test Tube Rack for 32 Test Tubes (16mm) is highly durable, and features adjustable angle and tube height support without requiring any base plate. It can be easily mounted and released.

Suitable for the following heating bath circulators: CW3-20, CW3-30, RW3-2025P, and RW3-3025P.



Specifications:

Diameter (mm)	16
Insert Capacity of Racks	CW3-20 CW3-30 RW3-2025P RW3-2035P RW3-3025P RW3-3035P
	2ea

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com