

## Jeio Tech Test Tube Rack for 50 Test Tubes (8mm)

**Product Code:** AAA41531

**Brand:** Jeio Tech

Jeio Tech's Test Tube Rack for 50 Test Tubes (8mm) is highly durable and feature adjustable angle and tube height support without requiring any base plate. They can easily be mounted and released.

Suitable for Jeio Tech 3.5 Litre Water Bath (BW-05H), Jeio Tech 3.5 Litre Water Bath (BW-05B), Jeio Tech 10 Litre Refrigerating & Heating Bath Circulator (RW3-1025P), and Jeio Tech 10 Litre Refrigerating & Heating Bath Circulator (RW3-1035P).



### Specifications:

<b>Diameter (mm)</b>	8
<b>Insert Capacity of Racks</b>	BW-05H BW-05B RW3-1025P RW3-1035P
	1 ea

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)