

## Jeio Tech Test Tube Rack for 86 Test Tubes (10mm)

**Product Code:** AAA44582

**Brand:** Jeio Tech

Jeio Tech's Test Tube Rack for 86 Test Tubes (10mm) is highly durable and feature adjustable angle and tube height support without requiring any base plate.

Suitable for Jeio Tech Heating Bath Circulators (CW3 Series) and Jeio Tech Refrigerating & Heating Bath Circulators (RW3 Series).

### Specifications:

<b>Diameter (mm)</b>	10
<b>Insert Capacity of Racks</b>	CW3-20 CW3-30 RW3-2025P RW3-3025P RW3-3040P
	2ea



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)