

Liebherr 360 Litre Spark-Free Laboratory Refrigerator (LKexv 3910 MediLine)

Product Code: LIEBLkexv3910

Brand: Liebherr

The Liebherr 360 Litre Spark-Free Laboratory Refrigerator (LKexv 3910 MediLine) features a digital temperature display that indicates the refrigerator's interior temperature to the exact degree, which can easily be read from the outside to give users immediate information and without the need to open the refrigerator.

Its alarm test function can be used to check whether an internally or externally connected alarm system is functioning correctly (refrigerator's performance is not affected).

Features

- Tested according to ATEX Directive
- Electronic control system with 0.1°C temperature adjustment
- Audible and visual alarms for door-open and temperature (user-adjustable)
- RS485 interface and volt-free contact for central monitoring
- Pre-drilled access ports



Specifications

Capacity (L)		360
Temperature Range (°C)		+ 3 to + 8
Air Circulation		Ventilated
Display		Digital
Defrost		Automatic
Alarm		Visual and audible (temperature and door)
Shelves		5 adjustable glass shelves
Doors		Self-closing solid door, right hand reversible hinge
Door Lock?		Yes
Material	Exterior	White Powder-Coated
	Interior	Moulded Polystyrol
Dimensions (W x D x H)	Internal	440 x 435 x 1635mm
	External	600 x 615 x 1840mm
Weight (kg)		68

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com