

Lec 130 Litre Basic Laboratory Refrigerator (LR111)

Product Code: LECLR111

Brand: Lec Medical

The Lec 130 Litre Basic Laboratory Refrigerator (LR111) has the power to easily slide under a general laboratory bench. Supplied with a door lock and two keys, this freezer provides laboratories with great security allowing only approved users to open it. Finished in white and spark-free, it has adjustable temperature control and auto defrost. It also contains three glass shelves and a useful drawer.



Features

- **Spark-Free Interior:** Offers safe and reliable storage of volatile liquids, reagents, chemicals, poisons, and samples with a spark-free interior that reduces the risk of internal explosions
- **Lock:** Security is essential in laboratories, so this freezer comes complete with a door lock and two keys, ensuring that only authorised individuals are able to access it
- Adjustable thermostat
- Auto defrost

Specifications

Net Capacity (L)	130
Temperature Range (°C)	2 to 10
Temperature Display	No
External Dimensions (W x D x H)	545 x 554 x 845mm
Lock	Yes
Shelves/Drawers	1/3
Defrost Capability	Auto Defrost
Spark-Free	Yes

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com