

Deltalab Microtube Racks with Lift-Off Telescopic Lids (for 1.5 or 2ml Tubes)

Product Code: M-410

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's Microtube Racks with Lift-Off Telescopic Lids (for 1.5ml or 2ml Tubes) are designed for long-term storage of samples at temperatures down to -90°C .

These laboratory racks can accommodate 50 microtubes of 1.5 or 2ml (hole diameter 12mm). They also feature four frosted side panels and a frosted rim on the edges of the lid provide additional labeling options.

The lids have a moulded-in alphanumeric location identification system with a clear centre panel for easy tube viewing.

Specifications:

Dimensions (l x w x h)	140 x 73 x 37.5mm
Material	Polypropylene
Sterility	Non-Sterile
Quantity	5



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com