

## Deltalab Microtube Racks with Lift-Off Telescopic Lids (for 0.5ml Tubes)

**Product Code:** M-450

**Brand:** Deltalab

Made of autoclavable polypropylene, Deltalab's Microtube Racks with Lift-Off Telescopic Lids (for 0.5ml Tubes) are designed for long-term storage of samples at temperatures down to  $-90^{\circ}\text{C}$ .

These laboratory racks can accommodate 50 microtubes of 0.5ml (hole diameter 10mm). They also feature four frosted side panels and a frosted rim on the edges of the lid provide additional labeling options.

The lids have a moulded-in alphanumeric location identification system with a clear centre panel for easy tube viewing.

### Specifications:

<b>Dimensions (l x w x h)</b>	120 x 60 x 37.5mm
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	5



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)