

## Deltalab Natural Polypropylene Racks with Lids (50 Holes)

**Product Code:** M-500

**Brand:** Deltalab

Deltalab's Natural Polypropylene Racks with Lids (50 Holes) are designed for storing samples down to 90°C. Autoclavable, their lids are jointed to the rack by hinges.

All of these laboratory racks are alphanumerically matrixed and designed for 1.5 and 2ml microtubes or 2.0ml cryotubes.



### Specifications:

<b>Colour</b>	Natural
<b>Holes</b>	50
<b>Dimensions (mm)</b>	71 x 153 x 53
<b>Material</b>	Polypropylene
<b>Pack Quantity</b>	20
<b>Sterility</b>	Non-Sterile

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)