

Deltalab Green PCR Workstation

Product Code: M-540

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's Green PCR Workstation is designed to work with different PCR volumes at the same time. It can accommodate up to 32 individual 0.2ml PCR tubes, 24 1.5/2ml and 16 0.5ml microtubes. Tube wells are alphanumerically matrixed with raised bold type for accurate sample identification. The rack also has a transparent lid for easy viewing of samples. Frosted writing panels on four sides allow labelling or writing. The rack is supplied with a removable rack accommodating 0.2ml PCR tubes.

Both racks withstand from -90°C to 121°C .

Specifications:

Dimensions	Main Rack	215 x 118 x 50mm
	Removable Rack	125 x 85 x 15mm
Colour	Green	
Material	Polypropylene	
Case Quantity	5	

