

## Deltalab Green PCR Workstation

**Product Code:** M-540

**Brand:** Deltalab

Made of autoclavable polypropylene, Deltalab's Green PCR Workstation is designed to work with different PCR volumes at the same time. It can accommodate up to 32 individual 0.2ml PCR tubes, 24 1.5/2ml and 16 0.5ml microtubes. Tube wells are alphanumerically matrixed with raised bold type for accurate sample identification. The rack also has a transparent lid for easy viewing of samples. Frosted writing panels on four sides allow labelling or writing. The rack is supplied with a removable rack accommodating 0.2ml PCR tubes.

Both racks withstand from  $-90^{\circ}\text{C}$  to  $121^{\circ}\text{C}$ .

### Specifications:

<b>Dimensions</b>	<b>Main Rack</b>	215 x 118 x 50mm
	<b>Removable Rack</b>	125 x 85 x 15mm
<b>Colour</b>	Green	
<b>Material</b>	Polypropylene	
<b>Quantity</b>	5	

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)