

Deltalab Yellow PCR Workstation

Product Code: M-542

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's Yellow PCR Workstation is designed to work with different PCR volumes at the same time. It can accommodate up to 32 individual 0.2ml PCR tubes, 24 1.5/2ml and 16 0.5ml microtubes. Tube wells are alphanumerically matrixed with raised bold type for accurate sample identification. The rack also has a transparent lid for easy viewing of samples. Frosted writing panels on four sides allow labelling or writing. The rack is supplied with a removable rack accommodating 0.2ml PCR tubes.

Both racks withstand from -90°C to 121°C .

Specifications:

Dimensions (mm)	Main Rack	215 x 118 x 50mm
	Removable Rack	125 x 85 x 15mm
Colour	Yellow	
Material	Polypropylene	
Quantity	5	

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com