

## Deltalab Pre-Assembled Freezing Cardboard Storage Box (M-600)

**Product Code:** M-600

**Brand:** Deltalab

Designed for the storage and preservation of frozen samples, Deltalab's Pre-Assembled Freezing Cardboard Storage Box (M-600) is made of white treated cardboard and covered by special white paper hand-stuck and hand-assembled for a better quality.

This laboratory storage box features pre-assembled compartments and can resist temperatures up to  $-100^{\circ}\text{C}$ . It is ideal for microtubes, cryovials, and the racks of freezer chests. It also has printed alphanumeric identifications and is silk screen printed on the base and lid.

### Specifications:

<b>External Dimensions (mm)</b>	150 x 55 x 55
<b>Numeric Grid</b>	10 x 10
<b>Holes (<math>\varnothing</math> mm)</b>	12
<b>Material</b>	Cardboard
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)