

Deltalab Pre-Assembled Freezing Cardboard Storage Box (M-600)

Product Code: M-600

Brand: Deltalab

Designed for the storage and preservation of frozen samples, Deltalab's Pre-Assembled Freezing Cardboard Storage Box (M-600) is made of white treated cardboard and covered by special white paper hand-stuck and hand-assembled for a better quality.

This laboratory storage box features pre-assembled compartments and can resist temperatures up to -100°C . It is ideal for microtubes, cryovials, and the racks of freezer chests. It also has printed alphanumeric identifications and is silk screen printed on the base and lid.

Specifications:

External Dimensions (mm)	150 x 55 x 55
Numeric Grid	10 x 10
Holes (Ø mm)	12
Material	Cardboard
Sterility	Non-Sterile
Quantity	1



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com