

Deltalab Freezing Cupboard Storage Box (134 x 134 x 57mm)

Product Code: M-606

Brand: Deltalab

Designed for storing and preserving frozen samples, Deltalab's Freezing Cupboard Storage Box (134 x 134 x 57mm) is made of white treated cardboard and covered by special white paper hand-stuck and hand-assembled for a better quality.

It has pre-assembled compartments and can resist temperatures up to -100°C .

Ideal for microtubes and cryovials, this storage box has printed alphanumeric identifications and is suitable for the racks of freezer chests.



Specifications:

| | |
|----------------------------|------------------|
| External Dimensions | 134 x 134 x 57mm |
| Grid Num. | 10 x 10 |
| Ø Holes (mm) | 10 |
| Material | Cardboard |
| Case Quantity | 10 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com