

Medline 16,000rpm Benchtop High Speed Centrifuge (MC-I6000)

Product Code: MC-I6000

Brand: Medline Scientific

The Medline 16,000rpm Benchtop High Speed Centrifuge (MC-I6000) is ideal for laboratory applications that need major flexibility, momentum, and processing. With its extensive range of motors and a maximum speed of 20,000rpm, this centrifuge provides a powerful solution to highly technical applications (e.g. specimen preparation) whilst taking up minimal bench space.



Specifications:

Type	High-Speed
Max. Speed (rpm)	16,000
Max. RCF (g)	17,800
Max. Capacity	12 x 10ml
Timer (min.)	0 - 99 min.
Motor	Brushless DC
Noise (dB)	≤ 55
Dimensions (l x w x h)	330 x 390 x 320mm
Weight (kg)	26

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com