

## Medline 20,000rpm Benchtop High Speed Centrifuge (MC-20000)

**Product Code:** MC-20000

**Brand:** Medline Scientific

The Medline 20,000rpm Benchtop High Speed Centrifuge (MC-20000) is ideal for laboratory applications that need major flexibility, momentum, and processing. With its extensive range of motors and a maximum speed of 20,000rpm, this centrifuge provides a powerful solution to highly technical applications (e.g. specimen preparation) whilst taking up minimal bench space.



### Specifications:

Type	High-Speed
Max. Speed (rpm)	20,000
Max. RCF (g)	27,800
Max. Capacity	4 x 100ml
Timer (min.)	0 - 99 min.
Motor	Brushless DC
Noise (dB)	≤ 55
Dimensions (l x w x h)	330 x 390 x 320mm
Weight (kg)	26

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)