

Medline 6000rpm Benchtop Low Speed Centrifuge (MC-7000)

Product Code: MC-7000

Brand: Medline Scientific

The highly adaptable Medline 6000rpm Benchtop Low Speed Centrifuge (MC-7000) has a variety of uses from routine specimen preparation to low-speed routine centrifugation, and is referred to by some manufacturers as universal due to its extensive range of uses and motors.

With its compact and robust design, it is simple to reach into the centrifuge whilst still taking up minimal bench space.



Specifications:

Type	Low Speed
Max. Speed (rpm)	6000
Max. RCF (g)	5080
Max. Capacity	4 x 500ml
Timer (min.)	0 - 99
Motor	Brushless
Noise (dB)	≤ 55
Dimensions (l x w x h)	440 x 550 x 330mm
Weight (kg)	45

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com