

Medical Disposable Protective Face Shield with Anti-Fog shield

Comfortable, latex-free full-face shield with elastic strap and foam forehead cushion.

Features:

- Comfortable, latex-free full-face shield with elastic strap and
- Double side anti-fog, anti-static, no-glare lens with the size of 32x24cm(12.6"x9.45").
- Hypoallergenic foam band can absorb sweat and provide enough room for eye glasses or safety goggles.
- Vented foam design for increased air flow and comfort.
- Sonically welded band gives the face shield added strength and reliability.
- Available in both elastic band and adjustable Velcro band.
- Lightweight and comfortable to wear, quick and easy to don.
- Excellent protection against potential contamination from blood pathogens, body fluids, or harmful chemical splash.



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com