

Medline 7000rpm Ultra Microcentrifuge (MDI008)

Product Code: MDI008

Brand: Medline Scientific

Medline's 7000rpm Ultra Microcentrifuge (MDI008) brings a new level of performance and reliability to laboratories. This mini microcentrifuge takes up very little space on the bench and is suitable for many applications including tissue/cell culture, molecular biology, clinical applications as well as proteomics and genomics. The DC motor produces a noise level of just 45 dB(A). Rotors are included.

Features

- Ideal for microfiltration and quick spin-downs
- Easy-to-use with dual start/stop function
- Fixed speed operation
- Quiet running
- Rotor change without tools
- Accelerates and brakes in seconds
- cTÜVus and FCC MCA certified
- Includes rotors

Specifications

Speed Range (rpm)	7000
Max. RCF (kg)	2680
Rotor Capacity	0.2/0.5/1.5/2mL x 8 0.2mL x 16 PCR strips or 0.2mL x 2 PCR 8 strips
Dimensions (D x W x H)	160 x 170 x 122mm
Weight (kg)	0.5



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com