

## Medline Mini Dry Block Heater (with heating lid)

**Product Code:** MD110401006

**Brand:** Medline Scientific

With its digital control wheel and bright LED display, Medline's Mini Dry Block Heater (with heating lid) is perfect for carrying out the incubation and activation of cultures, melting and boiling point determination, enzyme reactions, immunoassays, blood urea nitrogen determinations, and other laboratory applications.

Various interchangeable blocks and fixed microtube blocks available.

### Features

- Small footprint and compact structure
- Simple LED display for easy reading
- Heating lid provided
- Temperature is up to 100°C
- Simultaneous display of set/actual time and temperature
- 3D digital wheel adjustment of speed and temperature
- High temperature precision is up to  $\pm 0.5^{\circ}\text{C}$
- Programmable operation can store up to 5 programmes
- Low operation voltage (12V) can be used directly with car charger
- Various applications with interchangeable blocks or with fixed microtube blocks



## Specifications

<b>Temperature Control Range (°C)</b>	Ambient +5 to 100
<b>Timer</b>	0 - 99h. 59 min.
<b>Temperature Control Accuracy (°C)</b>	± 0.5
<b>Display Accuracy (°C)</b>	± 0.1
<b>Dimensions (W x D x H)</b>	164 x 116 x 98mm
<b>Weight (kg)</b>	0.9

## Accessories and Options:

Product Code	Description	Max. Mountable Capacity of Blocks
MD1102A	40 x 0.2ml Block	1
MD1102B	24 x 0.5ml Block	
MD1102C	15 x 15ml Block	
MD1102D	15 x 2ml Block	
MD1102E	8 - 12.5 x 12.5 x 32mm Block (used in cuvette)	
MD1102G	4 x 15ml Block	
MD1102H	2 x 50ml Block	

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)