

Medline 14 Litre Water Bath

Product Code: MDBS-14

Brand: Medline Scientific

Manufactured from polycarbonate with a safety thermostat, the Medline 14 Litre Water Bath provides optimal temperature control for samples.



Specifications:

Bath Volume (litres)		14
Temperature	Range (°C)	Amb. +5 to 100
	Accuracy (°C)	± 0.3 (at 37°C)
Dimensions	Bath Opening/Depth (w x l, d)	210 x 350, 160mm
	External (w x d x h)	280 x 500 x 300mm
Weight (kg)		6.5