

Medline 30 Litre Water Bath

Product Code: MDBS-30

Brand: Medline Scientific

Manufactured from polycarbonate with a safety thermostat, the Medline 14 Litre Water Bath provides optimal temperature control for samples.



Specifications:

Bath Volume (litres)		30
Temperature	Range (°C)	Amb. +5 to 100
	Accuracy (°C)	± 0.3 (at 37°C)
Dimensions	Bath Opening/Depth (w x l, d)	320 x 510, 180mm
	External (w x d x h)	380 x 660 x 240mm
Weight (kg)		11.2

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com