

## Medline 7 Litre Water Bath

**Product Code:** MDBS-7

**Brand:** Medline Scientific

Manufactured from polycarbonate with a safety thermostat, the Medline 7 Litre Water Bath provides optimal temperature control for samples.



## Specifications:

<b>Bath Volume (litres)</b>		7
<b>Temperature</b>	<b>Range (°C)</b>	Amb. +5 to 100
	<b>Accuracy (°C)</b>	± 0.3 (at 37°C)
<b>Dimensions</b>	<b>Bath Opening/Depth (w x l, d)</b>	160 x 300, 110mm
	<b>External (w x d x h)</b>	220 x 440 x 240mm
<b>Weight (kg)</b>		3.5