

Medline Polypropylene Universal Containers with Printed Label and Flow Seal Cap (30ml)

Product Code: 409205FL


Brand: Medline Scientific

Manufactured with ultra clear virgin polypropylene, Medline's 30ml Universal Containers can be used for collecting and storing a wide range of samples.

These containers are supplied sterile (inner bags of 100) and feature a printed label along with a leak-proof, flow-seal cap ensuring excellent sample containment.

Specifications:

Capacity (ml)	30
Material(s)	Polypropylene
Cap Colour	White
Sterility	Sterile
Quantity	400



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com