

Medline 0 - 1000g Precision Balance (TDI0002M)

Product Code: TDI0002M

Brand: Medline Scientific

With full capacity subtraction and overload protection, the Medline 0 - 1000g Precision Balance (TDI0002M) is designed for auto-calibration, routine weighing, and measuring density.

Equipped for laboratory applications in various sectors (e.g. industrial, education), this balance is exceptionally stable under all conditions and delivers accurate results.

Features:

- Auto calibration, check weighing, parts counting.
- High precision load cell.
- Full plastic housing, stainless steel pan.
- Super bright LCD display.
- Full capacity subtraction, overload protection.
- Below balance hanger.
- Stabilisation time is typically 1.5 seconds.
- Power supply, 9V DC/Dry battery.



Specifications:

Capacity (g)	0 to 1000
Readability (g)	0.01
Pan Size (mm)	Ø 133
Weighing Units	g, ct, oz
Weight (kg)	1.3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com