

## Medline 0 - 3100g Precision Balance (TD3000IM)

**Product Code:** TD3000IM

**Brand:** Medline Scientific

With full capacity subtraction and overload protection, the Medline 0 - 3100g Precision Balance (TD3000IM) is designed for auto-calibration, routine weighing, and measuring density.

Equipped for laboratory applications in various sectors (e.g. industrial, education), this balance is exceptionally stable under all conditions and delivers accurate results.

### Features:

- Auto calibration, check weighing, parts counting.
- High precision load cell.
- Full plastic housing, stainless steel pan.
- Super bright LCD display.
- Full capacity subtraction, overload protection.
- Below balance hanger.
- Stabilisation time is typically 1.5 seconds.
- Power supply, 9V DC/Dry battery.



## Specifications:

<b>Capacity (g)</b>	0 to 3100
<b>Readability (g)</b>	0.1
<b>Pan Size (mm)</b>	155 x 143
<b>Weighing Units</b>	g, ct, oz
<b>Weight (kg)</b>	1.3

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)