

Medline 0 - 5100g Precision Balance (TD5000IM)

Product Code: TD5000IM

Brand: Medline Scientific

With full capacity subtraction and overload protection, the Medline 0 - 5100g Precision Balance (TD5000IM) is designed for auto-calibration, routine weighing, and measuring density.

Equipped for laboratory applications in various sectors (e.g. industrial, education), this balance is exceptionally stable under all conditions and delivers accurate results.

Features:

- Auto calibration, check weighing, parts counting.
- High precision load cell.
- Full plastic housing, stainless steel pan.
- Super bright LCD display.
- Full capacity subtraction, overload protection.
- Below balance hanger.
- Stabilisation time is typically 1.5 seconds.
- Power supply, 9V DC/Dry battery.



Specifications:

Capacity (g)	0 to 5100
Readability (g)	0.1
Pan Size (mm)	155 x 143
Weighing Units	g, ct, oz
Weight (kg)	1.3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com