

## Medline 23 Litre Standard Autoclave

**Product Code:** MDB23C

**Brand:** Medline Scientific

The Medline 23 Litre Standard Autoclave is ideal for the modern laboratory in the pharmaceutical and food industries as well as for special applications. It comes with advanced features for safe, accurate, reproducible and validation for sterilisation in research and industry within the dental, medical, cosmetic, laboratory, veterinary, and tattoo fields. It is produced under the most rigorous system of quality control certifications in accordance to ISO13485 and TUV CE0197.

This lab autoclave is fully automated with an LCD display, an independent steam generator, temperature and pressure monitoring systems. This unit also comes equipped with a unique high precision control system and automatic alarm for water level, door locking and cycle errors, plus a built-in printer.



## Specifications:

<b>Volume (litres)</b>	23
<b>Chamber Size (mm)</b>	247 x 450
<b>Power (W)</b>	2200
<b>Printer</b>	Optional
<b>Double Door Locking</b>	Yes
<b>No. of Trays (max.)</b>	4
<b>Dimensions (w x d x h)</b>	480 x 455 x 700mm
<b>Weight (kg)</b>	58

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)