

## Medline 4.2kg Microincinerator

**Product Code:** BACTIPRO-500

**Brand:** Medline Scientific

The Medline 4.2kg Microincinerator is suitable for slide fixation, the sterilisation of loops, and culture tubes.

Its environmentally friendly design ensures minimum power consumption.

The operational angle adjustment of the chamber is 270°.



### Specifications:

<b>Energy Consumption</b>	0.150 kw/h
<b>Quartz Glass Tube Dimensions</b>	Ø 18mm, Length: 100mm
<b>Sterilisation Temperature</b>	850°C to 1100°C
<b>Time to Reach Sterilisation Temperature</b>	7 - 9 minutes
<b>Sterilisation Duration</b>	2 - 4 seconds
<b>Dimensions (w x d x h)</b>	200 x 300 x 265mm
<b>Weight (kg)</b>	4.2

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)