

Medline 8 Litre Standard Autoclave

Product Code: MDB8C

Brand: Medline Scientific

The Medline 8 Litre Standard Autoclave is ideal for the modern laboratory in the pharmaceutical and food industries as well as for special applications. It comes with advanced features for safe, accurate, reproducible and validation for sterilisation in research and industry within the dental, medical, cosmetic, laboratory, veterinary, and tattoo fields. It is produced under the most rigorous system of quality control certifications in accordance to ISO13485 and TUV CE0197.

This lab autoclave is fully automated with an LCD display, an independent steam generator, temperature and pressure monitoring systems. This unit also comes equipped with a unique high precision control system and automatic alarm for water level, door locking and cycle errors, plus a built-in printer.



Specifications:

Volume (litres)	8
Chamber Size (mm)	170 x 320
Power (W)	1500
Printer	Optional
Double Door Locking	Yes
No. of Trays (max.)	2
Dimensions (w x d x h)	414 x 365 x 570mm
Weight (kg)	34.5

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com