

Medline Bacti-Cinerator Sterilizer

Product Code: BACTI-100

Brand: Medline Scientific

Medline's Bacti-Cinerator Sterilizer sterilizes needles, loops, and culture tubes in 5 to 7 seconds at the optimum sterilizing temperature of 815.6°C (1500°F).

It has a ceramic funnel tube enclosed in a stainless steel (perforated) guard and a stainless steel support stand.



Specifications:

Dimensions (w x d x h)	130 x 102 x 180mm
Max. Diameter of Sterilizer Goods (mm)	14
Length (mm)	110
Weight (kg)	1.3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com