

Medline Blue Plastic Inoculating Loops (10 μ l)

Product Code: 302724R


Brand: Medline Scientific

Very flexible and sterile, Medline's Blue Plastic Inoculating Loops (10 μ l) provide regular wetting and full liquid transfer. Available with rounded tips, each inoculating loop has smooth edges and are all colour coded.

They are not recommended for colony counting.

Specifications:

Colour	Blue
Material	Polystyrene
Sterility	Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com