

## Medline Blue Sterile Smart Loops (5µl)

**Product Code:** 302715

**Brand:** Medline Scientific

Medline's Blue Sterile Smart Loops (5µl) are manufactured in flexible polystyrene and are sterilised by ethylene oxide.

Smooth and flexible, they facilitate uniform and smooth streaking without damaging the gel surface.


They are available in convenient sterile peel pouch bags with their batch and expiry date clearly displayed.

### Features:

- Rounded tip.
- Colour coded.
- Sterile and disposable to eliminate the risk of cross contamination.
- Free of any lubricants and do not contaminate any media.

### Specifications:

<b>Length (mm)</b>	173
<b>Weight (g)</b>	17
<b>Colour</b>	Blue
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)