

Medline Clear Sterile Smart Loops (1µl)

Product Code: 302713

Brand: Medline Scientific

Medline's Clear Sterile Smart Loops (1µl) are manufactured in flexible polystyrene and are sterilised by ethylene oxide.

Smooth and flexible, they facilitate uniform and smooth streaking without damaging the gel surface.

They are available in convenient sterile peel pouch bags with their batch and expiry date clearly displayed.

Features:

- Rounded tip.
- Colour coded.
- Sterile and disposable to eliminate the risk of cross contamination.
- Free of any lubricants and do not contaminate any media.

Specifications:

Length (mm)	173
Weight (g)	23
Colour	Clear
Material	Polystyrene
Sterility	Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com