

Medline L Shaped Spreaders (Pack of 25)


Product Code: MLSP25

Brand: Medline Scientific

Medline's L Shaped Spreaders (Pack of 25) are designed for surface spreading.

Their base length allows spreading to be carried out over any type of petri dish and avoids contact with the dish walls.

Available in convenient sterile peel-pouch bags, these spreaders are sterilized with ethylene oxide and their curved extremity helps to stop damaging the media surface while spreading.



Specifications:

Material	Polypropylene
Sterility	Sterile
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com