

## Medline Large PET Sweet Jars with Biohazard Labels (3600ml)

**Product Code:** 870036H

**Brand:** Medline Scientific

Used for the routine disposal of bench top waste (e.g. loops, cotton wool, wipes) and other small hazardous waste in every day microbiology or pathology procedures, Medline's Large PET Sweet Jars with Biohazard Labels (3600ml) are transparent and an excellent alternative to glass jars.

Biohazard labels are clearly displayed on the front of the jars.

### Specifications:

<b>Volume (ml)</b>	3600
<b>Height (mm)</b>	260
<b>Diameter (mm)</b>	140
<b>Neck Diameter (mm)</b>	140
<b>Material</b>	Polyethylene Terephthalate
<b>Quantity</b>	16

