

Medline Metal Cap Containers with Plain Labels (250ml)

Product Code: 309803M

Brand: Medline Scientific

Medline's Metal Cap Containers with Plain Labels (250ml) are sterile containers with a metal flowed seal inert liner cap.

Tray wrapped, these effective sample containers are manufactured from medical grade virgin polystyrene and the caps guarantee samples to be safely secured with their leak proof seal.

Specifications:

Volume (ml)	250
Labels	Plain
Material	Polystyrene
Sterility	Sterile
Quantity	50



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com