

## Medline Polypropylene Universal Containers with Printed Labels (30ml)

**Product Code:** 409201FL

**Brand:** Medline Scientific

Manufactured with ultra clear virgin polypropylene, Medline's Polypropylene Universal Containers with Printed Labels (30ml) can be used for collecting and storing a wide range of samples.

These containers are supplied sterile and feature printed labels, along with a leak-proof flow seal cap.

### Specifications:

Capacity (ml)	30
Material(s)	Polypropylene
Cap Colour	White
Sterility	Sterile
Quantity	400



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)