

Medline Polypropylene Universal Containers with Printed Labels & Boric Acid (30ml)

Product Code: 409202BP

Brand: Medline Scientific

Suitable for hospitals and clinics, Medline's Polypropylene Universal Containers with Printed & Boric Acid (30ml) contain a boric acid preservative that allows the storage and transport of urine samples to laboratories.

Their secure screw caps allow samples to be collected without causing any leakage.

Specifications:

Capacity (ml)	30
Material(s)	Polypropylene
Cap Colour	White
Sterility	Sterile
Quantity	400

