

Medline Polystyrene Universal Containers with Printed Labels & Boric Acid (30ml)

Product Code: 409202B

Brand: Medline Scientific

Supplied sterile, Medline's Polystyrene Universal Containers with Printed Labels & Boric Acid (30ml) contain a boric acid preservative that allows the storage and transport of urine samples to laboratories. All containers are manufactured from virgin, medical grade polystyrene.

Specifications:

Capacity (ml)	30
Material(s)	Polystyrene
Cap Colour	Red
Sterility	Sterile
Quantity	400



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com