

Medline Precision Balance (0 - 100g)

Product Code: MD8103

Brand: Medline Scientific

With a strong glass windshield and super bright LCD display with a back-light, the Medline Precision Balance (0 - 100g) is an excellent tool for carrying out general purpose weighing applications with its easy-to-use functions.

Features:

- High precision load cell.
- Aluminium alloy die cast base.
- Stainless steel pan.
- Super bright LCD display with backlight.
- AC/DC power supply.
- Full capacity subtraction.
- Overload protection.
- Glass shield supplied as standard.
- RS232 interface.
- Below balance hanger.
- Check weighing.
- Stabilisation time 1.5 seconds.



Specifications:

Capacity (g)	0 to 100
Readability (g)	0.001
Repeatability (g)	± 0.001
Linearity (g)	± 0.001
Pan Size (mm)	Ø 125
Weighing Units	g, ct, oz, lb, ozt
Dimensions (w x d x h)	190 x 287 x 267mm
Weight (kg)	1.3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com