

## Medline Precision Balance (0 - 200g)

**Product Code:** MD8104

**Brand:** Medline Scientific

With a strong glass windshield and super bright LCD display with a back-light, the Medline Precision Balance (0 - 200g) is an excellent tool for carrying out general purpose weighing applications with its easy-to-use functions.

### Features:

- High precision load cell.
- Aluminium alloy die cast base.
- Stainless steel pan.
- Super bright LCD display with backlight.
- AC/DC power supply.
- Full capacity subtraction.
- Overload protection.
- Glass shield supplied as standard.
- RS232 interface.
- Below balance hanger.
- Check weighing.
- Stabilisation time 1.5 seconds.



## Specifications:

<b>Capacity (g)</b>	0 to 200
<b>Readability (g)</b>	0.001
<b>Repeatability (g)</b>	± 0.001
<b>Linearity (g)</b>	± 0.001
<b>Pan Size (mm)</b>	Ø 125
<b>Weighing Units</b>	g, ct, oz, lb, ozt
<b>Dimensions (w x d x h)</b>	190 x 287 x 267mm
<b>Weight (kg)</b>	1.3

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)