

## Medline Rotary Vane Vacuum Pump (CCMD4)

**Product Code:** MED-141C

**Brand:** Medline Scientific

Equipped with a direct connect oil seal to prevent air reflux, Medline's Rotary Vane Vacuum Pump (CCMD4) is specially designed for vacuum ovens to achieve a higher level of vacuum pressure.

Ideal for low laboratory applications.

An optional oil mist trap is available to prevent oil in-flow. Recommended for use with Jeio Tech's Vacuum Ovens (OV Series).

### Specifications:

<b>Stages</b>	2
<b>Pumping Speed (L/min.)</b>	140
<b>Ultimate Vacuum (m/bar)</b>	5 x 10 <sup>-3</sup>
<b>Gas Ballast</b>	Yes
<b>Motor Power (W)</b>	370
<b>Max. Water Vapour Pumping Rate (kg/hr.)</b>	400
<b>Rotational Motor Speed (rpm)</b>	2880
<b>Oil Capacity (ml)</b>	450
<b>Recommended Oil (litres)</b>	1.1
<b>Inlet/Outlet Port (DN)</b>	8mm Hose Tail
<b>Oil Mist Eliminator (Exhaust Outlet)</b>	Fitted
<b>Power Compatibility (Single Phase)</b>	110/240V Dual
<b>Dimensions (w x l x h)</b>	178 x 375 x 330mm
<b>Weight (kg)</b>	12



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)