

## Medline Rubber Drain Tube (6.3 x 1.7mm)

**Product Code:** 211492

**Brand:** Medline Scientific

The Medline Rubber Drain Tube (6.3 x 1.7mm) is manufactured from 45° natural rubber, high strength, good resistance to tear, fatigue, and abrasion. Has excellent electrical properties but poor resistance to ozone, mineral, and petroleum based oils.



### Specifications:

|                        |  |
|------------------------|--|
| <b>Temperature</b>     | -65°C to +70°C                                 |
| <b>Applications</b>    | Transferring gases, liquids, and diluted acids |
| <b>Coil Length (m)</b> | 20   |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)